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
asatim

angled solid abutment treatable in mouth

TRAINING COURSE FOR THE ASATIM SYSTEM
**NEW SYSTEM OF ABUTMENTS FOR CEMENTED PROSTHESES SUPPORTED
BY DENTAL IMPLANTS**

Training course for the ASATIM system

TRAINING COURSE FOR THE ASATIM SYSTEM

Abutments for cemented prostheses supported by dental implants

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Organization and reservations:
610898503 – (Ask for Fortu)
963302461 – Main Office
tecnicogarante@asatim.net

Schedule: Friday 10 am -2 pm and 4 pm-8 pm. Saturday 10 am -2 pm
A maximum of 8 places are available per course.

At the completion of the course, participants will receive a diploma issued by **ASATIM S.L.**

Documents pertaining to the **ASATIM** system:

- European Patent EP-1772113
- US patent US20070082320A1

CE - 0476

ISO 13485

ISO 9001

SCOPE OF THE **ASATIM** COURSES:

- General information about **ASATIM** solid straight and angled abutments.
Using the asatim box.
- Practical handling of models in an adjustable articulator: selecting the appropriate **ASATIM** abutments.
- Sealing the **ASATIM** abutments on implants in the articulated model.
- Trimming the **ASATIM** abutments to achieve parallelism among them so they can support a fixed cemented prosthesis.
- Adjusting a pre-formed resin provisional on the **ASATIM** abutments.
- Rethreading an implant where the thread of a broken abutment is lodged inside.

Training course for the ASATIM system

THE ASATIM SYSTEM

The ASATIM system is a procedure for placing and managing one-piece solid straight and angled abutments on osseointegrated dental implants. The procedure is done entirely at the dental office, avoiding any need to use the dental laboratory for the process.

Most dentists agree that the surgical phase is often simpler than the prosthetic phase, primarily due to the many steps and tools involved in the latter.

The **ASATIM** system is a prosthetic method that will give you great autonomy and freedom in the methodology and manufacture of fixed dental prostheses cemented on implants.

Using **ASATIM** abutments, you will be able to place a provisional prosthesis very economically, until the definitive ceramic prosthesis is ready, thus reducing considerably the total cost of each prosthesis.

Worldwide, there are two main types of connections in the dental implant market: Compatible Internal Hexagon (CIH), and Compatible External Hexagon (CEH). These connections are made by approximately 60-70% of all implant companies.

The ASATIM system manufactures compatible abutments for Compatible Internal Hexagon (CIH, Niznick), and for Compatible External Hexagon (CEH).

In the procedures for oral rehabilitations on implants, the ideal situation is to place parallel implants in the jaw, so that the prosthetic abutments used on them can be straight, but this is not always possible.

There are different methods or systems to solve the problem of implant unparallelism when making a cemented fixed dental prosthesis. This is the main focus of the **ASATIM** system.

Example of inclined implants avoiding the maxillary sinus.

Depending on the type of dental prosthesis supported by the implants, there are several systems:

- prostheses which are threaded to the implant directly or with multi-units.
- removable prostheses supported by implants and the gingiva.
- fixed prosthesis cemented to already parallelised abutments.

The ASATIM system is designed for cemented prostheses.

Professional practitioners have two options to achieve parallelism of the abutments so they can support a cemented prosthesis:

- 1) Let the dental laboratory solve the unparallelism.

With this option, the dentist's tasks are limited to the surgical placement of the implants. For the prosthetic phase, the dentist simply makes impressions to transfer the working model to the dental laboratory, and leaves to the dental prosthetic technician the problem of overcoming the unparallelism.

- 2) Try to solve the unparallelism in the dental office. Using the abutments provided by his/her company of choice, the practitioner must correct implant unparallelism. However, when trying to parallelise the hollow abutments provided by the chosen company, which feature an anchoring screw, it becomes obvious that, especially when the number of abutments increases, some areas of the abutments cannot be trimmed or filed sufficiently, because they are hollow.

How does ASATIM overcome the problem of trimming hollow abutments?

The ASATIM system overcomes the problem by offering the possibility of trimming SOLID straight and angled abutments directly in the patient's mouth, using abundant irrigation, so that abutment unparallelism can be solved.

What is the ASATIM system?

The ASATIM system is a procedure for preparing, placing and managing straight and angled solid prosthetic abutments on osseointegrated dental implants. These abutments are processed entirely in the dental office, thus avoiding the necessary step of resorting to the dental laboratory.

Training course for the ASATIM system

The **ASATIM** system has been developed to prevent many of the causes of prosthetic failure, including:

- Inaccurate planning of the number and size of the implants.
- Failure to prevent seepage into the implant-abutment connection.
- Improper adjustment to the anatomy and physiology of the gingiva.
- Patient dissatisfaction due to aesthetic problems with the prosthesis.

What are the components of the Asatim system?

The **ASATIM** system comprises the following components:

- 1) The Asatim box
- 2) Straight and angled prosthetic abutments
- 3) Tightening keys
- 4) Instructions for selecting the appropriate angled abutments

1. The **ASATIM** Box:


It is a circular box (16 cm diameter) featuring a central axis or **connection hub** that simulates the connection of the implant.

The circular box is made of a sanitary bioplastic material. It features 12 compartments, distributed along the circumference of the box, and each compartment occupies 30° of the 360° space. Each compartment represents a **rotational position**, or angular orientation. Each one of the 12 compartments is divided in two parts to store and classify the different prosthetic abutments according to their angle of tilt, 14° or 24°. The box has a revolving lid that allows accessing the abutments.


The box can be sterilised.

The basic set of **ASATIM** abutments comprises:

- 12 solid abutments angled at 14°
- 12 solid abutments angled at 24°

Connection hub of the **ASATIM** box

Body of the **ASATIM** box with its 12 double compartments, without the lid

The figure shows the complete **ASATIM** box, with the rotating lid, and a few **ASATIM** system solid angled abutments, organised in their respective compartments

2. Angled prosthetic abutments:

The solid, one-piece angled or straight prosthetic abutments are designed to be trimmed directly in the mouth using abundant irrigation. The starting point of the thread in each **ASATIM** solid angled abutment determines the segment, within a circumferential space, where the tip of the angled portion of said abutment points when it is fully tightened into the central connection hub of the **ASATIM** box. This is the angled abutment's rotational position.

From left to right, Compatible Internal Hexagon (CIH) Abutment. Ø: 3.5 mm.

- CIH straight conical
- CIH solid angled at 14°
- CIH solid angled at 24°

From left to right, Compatible Internal Hexagon Abutment-Large (CIH-G). Ø: 3.9 mm.

- CIH-G, straight conical
- CIH-G solid angled at 14°
- CIH-G solid angled at 24°

From left to right, Compatible External Hexagon Abutment (CEH). Ø: 3.9 mm.

- CEH straight conical
- CEH solid angled at 14°
- CEH solid angled at 24°

3. Tightening keys:

Short hexagonal keys are provided to manage the abutments. These are the best-suited keys, due to the eccentric screwing action.

Solid abutments (CEH) angled at 24° inserted in the short hexagonal key (ES: 1.25)

4. Instructions for use:

A practical/instructional course is offered to learn how the **ASATIM** system works. The course takes place in adapted classrooms with dental models and their corresponding implants. Using techniques that help deduce the rotational position, participants learn to select the appropriate **ASATIM** solid angled abutments, to then screw them into implants placed in the dental arch of an articulated model.

Sealing and trimming is also practiced to achieve parallelisation, so that the abutments can support a cemented prosthesis.

Training course for the ASATIM system

Advantages and disadvantages of cemented restorations on “anchoring-screw abutments” vs. “screw-on restorations”

Advantages

- a) Cemented suprastructures are more passive due to the space they leave to allow for the cement. On the other hand, the anchoring screws of screwed prostheses never go in passively.
- b) The occlusal surfaces of cemented structures stay intact because they do not include orifices for the anchoring screws, allowing a more favourable axial load.
- c) There is a lower incidence of fracture of the occlusal porcelain in cemented prostheses due to the integrity of the occlusal surface.
- d) The manufacturing techniques used at the dental office are practically the same as those used to make conventional fixed prostheses.
- e) Manipulations in the posterior region are easier when cement is used, compared to screwed prostheses and using the screw driver. Cementing is a common, well-known dentistry process.
- f) Restorations using screwed prostheses are always more expensive, because more components are required, as well as more adjustments at the dental office and the laboratory.
- g) Supplying the small parts for screwed prostheses is more complex and costly.
- h) Prostheses have to be removed for check ups and maintenance. The follow-up for cemented prostheses focuses on checking the clinical appearance of the gingiva (for bleeding, halitosis, inflammation) and on making follow-up X-rays.
- i) If the porcelain fractures, it can be restored using traditional techniques, although, theoretically, with screw-on restorations it can be made at the dental laboratory.
- j) The aesthetic outcome of cemented prostheses is superior, since the heads of the screws, or the pads used to cover them, do not disfigure the porcelain surfaces. The figure shows such orifices in the structure at 22 and 23 to allow passage of the retaining screws from the structure to the implants.
- k) The figure also shows that the hexagons, which should fit into the implants, have disappeared, thus reducing the strength of the structure, because the retention depends solely on the anchoring screw.

Explanatory figure Advantage B

Explanatory figure Advantage G

Explanatory figure Advantage J

Explanatory figure Disadvantage C

Disadvantages

- a) A tight occlusal space forces the reorganisation of the opposite occlusal area, which usually is extruded and requires levelling.
- b) Special care is required with cement residue, which sometimes may require using anaesthesia to be removed.
- c) Both in cemented prostheses and in screw-on prostheses, **loosening of the anchoring screws increases the mobility of the components of the implants**, leading to fracture of the screws inside the implants, fracture of the porcelain, and deterioration of the internal threads of the implants.

However, this problem is easier to solve in screw-on restorations, because it is possible to screw the prosthesis again through the orifices in the crowns.

On the other hand, in cemented restorations, the prosthesis has to be removed by roughing it down with the air hand piece and abundant irrigation, which is a difficult task to perform.

This is the main problem that may occur, and it can happen rather frequently.

Perhaps because of this problem, dentists are leaning towards screw-on prostheses.

The ASATIM system is a newly available procedure designed to be used by itself in rehabilitations, or as a complement to cemented prosthetic systems already available in the market.

The different manufacturing companies offer straight or angled abutment pillars with an anchoring screw. This forces the abutments to be hollow, implying that they cannot be trimmed or roughed down to achieve parallelism among them, which very often would be desirable. However, using them in combination with the ASATIM system, their use options increase tremendously, because solid abutments can be paralleled to the hollow abutments.

Advantages and strong points of the ASATIM system vs. traditional screw-retained prostheses and cemented prostheses (hollow abutments with anchoring screw)

Advantages of the ASATIM system

- a) The ASATIM system is a technique that is done exclusively at the dental office, which saves time and requires fewer patient visits.
- b) It eliminates the time required by the dental laboratory.
- c) It eliminates transfer abutments to carry the position of the implants to the plaster model, including the transfer, the analogue, the laboratory anchoring screw, calcinables, soldering and taking silicon impressions.
- d) It prevents the inaccuracies which are typical of the process of transferring to the laboratory. The fitting of the prosthesis becomes simpler.
- e) High quality treatments are made while greatly reducing costs.
- f) Preparing the abutments inside the oral cavity is more accurate than doing it at the dental laboratory. Also, ASATIM abutments have thicker walls, allowing great freedom during their preparation.
- g) With the ASATIM system **the abutments never loosen or unscrew in their connection to the implant**, since the fixed cemented prosthesis keeps them from turning.

On the other hand, ASATIM abutments are not meant for single implant prostheses, since they do not have an anti-rotational system and, over time, the abutment can turn with the cemented crown, and the threaded portion of the ASATIM abutment can fracture.



STRONG POINTS:

- a) Sealing prevents filtration.
- b) Eccentric screwing is possible with a solid abutment.
- c) Intraoral trimming allows more accurate parallelism.

Training course for the ASATIM system

A) SEALING

In the **ASATIM** system, sealing is done to keep the abutment from unscrewing during the intraoral trimming, and to prevent seepage/filtration between the abutment and the implant.

These figures show the space of 85 microns, which is sufficient for the passage of biomolecules and bacteria.

Sealing is extremely important. It makes the difference between success and failure.

Most implant-abutment sets in the market are made with exceptional accuracy (if the calcinable system is not used it can be up to microns), but they operate in a septic environment, which implies the penetration of bacteria and lipid biomolecules from the first day they are placed, unless proper sealing is done.

This “greasing” process, together with the vibrations from the chewing cycles, but especially from bruxism and tooth tightening habits, lead to screw loosening, and ultimately to fracture, if proper sealing is not carried out.

ASATIM proposes using dual-cure resin (DCR). The fluid obtained after mixing the two components is introduced in the thread, in the internal or external geometric shape, in the metal-metal interface, and any other openings there may be, sealing everything to prevent any filtration/contamination.

When DCR is applied on an implant in the oral environment, the inside of the connection has to be dried first, and then the adhesive mixture is introduced using a fine brush. Immediately after, the prosthetic abutment is screwed on (if it is a one-piece abutment), or the anchoring screw with the hollow pillar is put in place.

Example of how to make the seal with dual cure resin. Wet the thread of the attachment with dual cure resin and also insert it with a brush inside the implant connection.

There is enough time to carry out this procedure before the adhesive begins to polymerise. If, for any reason, the process has to be interrupted after applying the DCR (the patient coughs or makes a sudden movement), the inside of the connection must be cleaned with the water spray and the operation can be repeated without any problem.

Interposing the dual-cure resin binds the implant to the solid prosthetic abutment, or to the hollow prosthetic abutment with the anchoring screw, so that it becomes one piece.

The best implant-prosthetic abutment connection would be a non-existing one (i.e., one-piece implants), then there would not be any risk of contamination/filtration. Introducing dual-cure resin is the most similar result we can achieve.

As we saw, interposing dual-cure resin binds the implant to the solid prosthetic abutment, or to the hollow prosthetic abutment with the anchoring screw, so that it becomes one piece, preventing their separation. As a result, one of the most frequently asked questions is how to solve the different complications that may arise in a prosthetic rehabilitation supported on implants.

What should be done when an intermediate implant fails in an implant- supported rehabilitation? (Figure A)

To solve this problem it has to be approached as if it were the surgical procedure for an apicoectomy: the flap is lifted with a pericoronal incision, and the granulomatous tissue that surrounds the implant is located; it is properly exposed, and then the implant is separated from the prosthesis using a knife-edge bur, just above the crown, being careful not to touch the ceramic; then the implant is removed, and it will come out easily; next the cyst cavity that usually forms around it is cleaned, it is filled with biomaterials if needed, and closed with a stitch.

0

1

2

3

FIGURE A | Figure A (0) shows the orthopantomography of an oral rehabilitation on implants. The failed implants to be removed are marked in red: a mesial implant in the mandible, and an intermediate one in the maxilla. Figure A (1) shows the mesial implant in the rehabilitation which will be removed. Figure A (2) shows the cut made in the lower prosthesis. Figure A (3) shows the removed implant

Training course for the ASATIM system

4

5

6

7

FIGURE A (continued) | Figure A (4) shows the crown (13) supported by the failed implant which will be removed from the oral rehabilitation in the maxilla. Figure A (5) shows the apical recession flap where the failed implant is visible. Figure A (6) shows the implant already cut at the gingival edge of the ceramic, being separated from the prosthesis. Figure A (7) shows the removed implant.

8

9

10

FIGURE A (continued) Figure A (8) shows the granuloma surrounding the failed implant. Figure A (9) shows the cavity which is left after the implant bed has been cleaned. Figure A (10) shows the suture.

What should be done when a mesial or distal implant fails in an oral rehabilitation?

If the failed implant is mesial or distal in a rehabilitation on implants, it should be separated from the prosthesis by cutting the prosthesis, so that the implant is freed, and then it can be separated easily.

What should be done when the ceramic fractures in an oral rehabilitation?

If the complication is the fracture of the ceramic of a crown, as it can happen with any other type of rehabilitation, our method is to fully trim the crown using the $\frac{3}{4}$ technique, without affecting the existing metal, and making a zirconium crown.

11

12

13

FIGURE A | Figure A (11) shows the fracture of the ceramic of a crown (23) that is part of an oral rehabilitation in the maxilla. Figure A (12) shows the trimming of the crown using the $\frac{3}{4}$ technique, and leaving the structure's metal. Figure A (13) the $\frac{3}{4}$ crown zirconium/ceramic already cemented.

B) ECCENTRIC SCREWING. TIGHTENING KEYS FOR THE ASATIM PROSTHETIC ABUTMENTS

Eccentric screwing: using the different hexagonal keys already available in the market (ES: 1.25). The key is inserted in the top end of the angled pillar. Since the screwing action is eccentric, the best suited keys are the short ones.

0

1

FIGURE B | Figure B (0) shows the **ASATIM** box, with a solid angled abutment screwed into the implant replica in the centre of the box. Its rotational position is 180°. The short hexagonal key is still inserted in the abutment (ES: 1.25). Figure B (1) shows a solid angled abutment, with a compatible internal hexagon connection, with the short hexagonal key inserted in it. The image also shows a solid angled abutment inside an implant with a compatible internal hexagon connection, making a switching-platform effect.


C) TRIMMING INSIDE THE MOUTH

After sealing the abutments, we proceed to trimming with abundant irrigation.

- Usually we start with the vertical dimension by cutting the height of the abutments, so that the patient feels comfortable as soon as possible and is able to close his/her mouth. The ideal height for the abutment, to provide optimal retention of the prosthetic crown and facilitate an esthetical prosthesis, is 5-7 mm from where it emerges in the gum. Knife-edge diamond burs are used for this procedure.
- Once the vertical dimension is cut to the desired height, the abutments are trimmed, or parallelised, so they can support a cemented prosthesis. Chamfer diamond and tungsten carbide burs are used for this procedure.

Training course for the ASATIM system

The advantage of the system comes from working with a solid abutment, which allows a significant removal of the titanium material without weakening the abutment or reducing its ability to retain the prosthetic crown (Figure C).

0

1

2

FIGURE C | Figure C (0) shows an orthopantomography where an oral rehabilitation will be done in the jaw with the lower front exodonty to place 6 implants, which will be very unparallel because of the lack of space in the bone. Implants will not be placed over the teeth; rather, the posterior limiting implant will be very tilted to distalize the emergence of the implant as distally as possible. This results in great unparallelism in the emerging abutments. Figure C (1) shows the orthopantomography of the performed treatment. Figure C (2) shows the emergence of the straight conical abutments, where their unparallelism is very noticeable. The parallelisation, or trimming, of the abutments is done as with natural teeth.

3

4

5

FIGURE C (continued) | Figure C (3) shows how the most appropriate **ASATIM** solid angled abutments have been selected. Next, the abutments are trimmed to be parallel among themselves to receive the prosthetic structure. Since they are solid, a lot of the material can be removed. Figure C (4) shows the achieved parallelism of the solid abutments. Figure C (5) shows the cemented metal-ceramic structure.

(The parallelisation or trimming of the abutments is managed as if they were normal teeth.)

Cases of rehabilitations involving 8-10 implants per dental arch (Figure D):

The goal is to parallelise the pillars, regardless of whether we have to rehabilitate a class II or class III occlusion, since the malocclusion compensation has to be done with the metal ceramic structure. Therefore, it is not necessary to trim a tilted shape to compensate a class III, even if it were parallel. It is better to parallelise perpendicularly to the gingiva (Figure D).

FIGURE D | Figure D (0) shows an orthopantomography where we can see 8 implants in the maxilla. They are very unparallel. Although it is not visible in this image, they are class III for the jaw rehabilitation. Figure D (1) shows the emergence of the solid angled abutments, already parallelised to each other to receive the prosthesis. Their emergence is perpendicular to the gingiva, with no direction compensations for class III. This means the compensation for the class III occlusion has to be compensated by the metal-ceramic structure itself. Figure D (2) shows the class III discrepancy compared to the lower rehabilitation. Figure D (3) shows the finished oral rehabilitation, in which the upper metal-ceramic structure is doing the class III occlusion compensation.

asatim

angled solid abutment treatable in mouth

SUMINISTROS ASATIM S.L.
C/ Mas del Bombo 4-A 46530 Puzol
Valencia, Tfno: 963302461

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***HEC: Conexión Europea del Dr Branemark (Plataforma Ø4,1) Hexágono externo compatible a las siguientes marcas comerciales registradas;** (European connection of Dr Branemark (Platform Ø4,1) External hexagon compatible with the following registered trademarks):
***HIC: Conexión americana del Dr Niznick (Plataforma Ø3,5) Hexágono interno compatible a las siguientes marcas comerciales registradas.** (American connection of Dr Niznick (Platform Ø3,5) Internal hexagon compatible with the following registered trademarks):